

# FORE-SIGHT® Roll Stand

Roll Stand with Basket and Padded Grip Handle



**FORE-SIGHT**  
Absolute Tissue Oximeter

**CAS Catalog Number: 01-06-0238**

**Specifications:**

- Load capacity: 30 lbs (13.6 kg)
- Overall height to mounting plate: 43" (109 cm)
- 21" (53.3 cm) diameter aluminum base
- 4" (10 cm) dual wheel casters (2 locking, 3 non-locking) for easy transition over uneven surfaces
- Built-in 10 lbs (4.5 kg) counterweight (beneath the base) for tip safety
- 5" (12.7 cm) Slide – In style mounting plate
- Adjustable tilt angle for optimum viewing
- Utility basket with hooks for cable management  
12" w x 6" h x 8" d (30.5 x 20.3 x 15.2 cm)
- Padded grip handle
- Requires FORE-SIGHT adapter plate  
P/N 01-06-0237

# FORE-SIGHT Mounting Accessories



**01-06-0237**  
Adapter Plate



**01-06-0223**  
FS 6" Utility Basket  
For attachment to 16" & 8" Arm wall mount  
and 16" Arm pole mount



**01-06-0232**  
FS 3" Utility Basket  
For attachment to 16" & 8" Arm wall mount  
and 16" Arm pole mount



**01-06-0236**  
FS 16" Arm Pole Mount  
Adapter Plate to be purchased separately



**01-06-0235**  
FS 8" Arm Wall Mount  
(Giraffe Bed)  
Adapter Plate to be purchased separately



**01-06-0234**  
FS 16" Arm Wall Mount  
Adapter Plate to be purchased separately



**01-06-0231**  
Giraffe Dovetail Attachment

**FORE-SIGHT**  
Absolute Tissue Oximeter

For ordering information, contact Customer Support: [custrv@casmed.com](mailto:custrv@casmed.com)  
CAS Medical Systems Inc. 44 E Industrial Rd, Branford, CT 06405 USA  
Toll Free: 800.227.4414 Phone: 203.488.6056

Visit us on the web at [www.cerebraloximetry.com](http://www.cerebraloximetry.com)