

PVC Free Blood Pressure Cuffs

Now available in two styles!



SoftCheck Free[®] Neonatal BP Cuffs

UltraCheck[®] BP Cuffs*

Why PVC Free?

- During the manufacture and incineration of PVC, Dioxin, a known carcinogen may be released into the atmosphere.
- DEHP, a phthalate used to soften PVC plastics can leach from medical devices and is linked to reproductive birth defects and other illnesses.



CASMED supplies two product lines devoted specifically to the protection of immuno-suppressed/ high-risk patients by offering PVC/DEHP Free blood pressure cuffs.

*UltraCheck manual cuffs are sold with a PVC, but DEHP-free handbulb.

CASMED[®]
FOR WHAT'S VITAL